

# LAr Filter Regeneration Procedure

Yun-Tse Tsai

April 17, 2022

Notes:

- P&ID is attached in the end of this list
- The weight for the all metal valve, V9, is 4 pounds
- The tube upstream V9 is not leak tight (mainly the connections to the flowmeter). Therefore keep V9 closed when we need to isolate the LAr filter
- The gas flow has to be greater than 130 scfm (marked on the flowmeter) to prevent the heater from getting too hot
- We should keep the gas flow (Ar or 2%H<sub>2</sub>+Ar) between 130 and 400 scfm (label on the flowmeter), but preferably at 400 scfm (the upper range of the flowmeter)
- Maintain the catalyst temperature between about 170°C and 225°C
- Do NOT exceed 225°C, even though 250°C may be tolerated
- If we keep the gas flow at 400 scfm, a gas bottle should be finished in about an hour
- If seeing smoke or smelling something unusual, shut down the variac power supply (heater) and investigate.

Checklist	What to Do and Detailed Description
<b>Preparation</b>	
3 bottle of ultra high purity Ar gas (TBC)	
5 bottles of Ar+2%H <sub>2</sub> gas (TBC)	
Tubes connecting the heater and the LAr filter wrapped with aluminum foils for thermal insulation	
V4 connected to the scroll pump	Prepare to evacuate the vacuum vessel
V4 open, scroll pump on	Evacuate the vacuum vessel
V3, V5, V6, V7, V8, V9, V10, V11, V12 closed	
Exhausting gas line connected and humidity meter hooked	
<b>Preheating with Ar gas</b>	
V3, V5, V6, V7, V8, V9, V10, V11, V12 closed	
PG6 at 0 psi	
Variac power supply off. Voltage set at 0	

Heater plugged in to the variac power supply	
Ar gas bottle connected to Reg1 and V7/V8 line	
GMV1 opened, Reg1 increased, V7 opened, air purged	Purge the air in the connection tube
V7 closed	Finish purging
V8, V9, V6 opened	
PG3 > 2 psig, V10 opened	
Gas flow > 130 scfh, stable	
Variac power supply on, increase the voltage	Turn on the heater
Humidity plateaued for > 10 minutes	Molecular sieves regenerated
Preheated for > 2 hours	
Variac power supply off. Voltage set at 0	Turn off the heater
V8, V9, V10 closed	
GMV1 and Reg1 closed	
<hr/>	
<b>Regenerating copper sieves</b>	
Ar+2%H <sub>2</sub> gas bottle connected to Reg1 and V7/V8 line	
GMV1 opened, Reg1 increased, V7 opened, air purged	Purge the air in the connection tube
V7 closed	Finish purging
V8, V9, V6 opened	
PG3 > 2 psig, V10 opened	
Gas flow between 50 and 160 slpm (Ar), or between 130 and 400 scfh (labeled as Air). Preferably at 400 scfh Air	
Variac power supply on, increase the voltage	Turn on the heater
Should the temperature exceed 225°C anywhere in the bed, switch to H <sub>2</sub> -free gas until the hot zone cools back down to 200 – 210°C, then resume feeding the H <sub>2</sub> gas mixture	
The temperature of the all catalyst bed is stable or subsiding	
Humidity plateaued for > 10 minutes	Copper sieves regenerated
Variac power supply off. Voltage set at 0	Turn off the heater
V6, V8, V9, V10 closed	
GMV1 and Reg1 closed	
Humidity meter unmounted	
<hr/>	
<b>Completion</b>	

V3, V5, V6, V9 closed

Scroll pump connected to V5

V5 open, scroll pump on

PG3 at 0 psig

V5 closed

V6, V10, V11, V12 closed

Prepare to evacuate the LAr filter

Scroll pump connected to V10

V10 open, scroll pump on

Pumping 10 minutes

V10 closed

V3, V5, V9 closed

V17 closed, turbo pump off?

Prepare to evacuate the piece of the plumbing system

V12, V6 open?

---




Figure 1: P&ID